



# Development of a human T cell backbone flow cytometry panel enabling flexibility in reagent choice while minimizing panel design challenges

Xiaoshan Shi,<sup>1</sup> Chad Sisouvanthong,<sup>2</sup> Bipulendu Jena,<sup>2</sup> Stephanie Widmann,<sup>2</sup> Aaron Tyznik<sup>2</sup>  
<sup>1</sup>BD Biosciences, San Jose, CA; <sup>2</sup>BD Biosciences, San Diego, CA

## Abstract

Expansion of existing flow cytometry panels with new markers of interest can result in suboptimal resolution and, in some cases, the need to design a new panel, impacting cost and increasing time to insight. To minimize these challenges and to provide increased flexibility, we have developed a human T cell backbone panel strategically designed to be complemented with 4–5 drop-in fluorochromes and markers of choice, depending on instrument configuration, with minimal panel design effort. The backbone panel contains five T cell markers (CD3, CD4, CD8, CD45RA, CCR7) conventionally used to identify different maturational states of CD4+ and CD8+ T cells (naïve, central memory, effector memory, effector memory RA). We will show how this backbone panel meets four fundamental requirements: i) clear resolution of major T cell subsets; ii) the fluorochromes used in the backbone panel have minimal resolution impact on the detectors allocated for drop-ins; iii) the fluorochromes assigned to drop-ins have minimal impact on the resolution of the backbone panel; iv) the fluorochromes assigned to drop-ins do not impact each other. We will provide examples of different drop-in combinations for the study of T cell and regulatory T cell biology. Performance of the backbone panel on different instruments will be shown to demonstrate assay consistency across platforms. The compatibility of the backbone panel with intracellular stain and transcription factor analysis will be further demonstrated.


Backbone		3–4 laser flow cytometer	5 laser flow cytometer
Marker	Fluorochrome	Drop-in fluorochromes	Drop-in fluorochromes
CD3	BV510	BV421/V450/	BV421/V450/
CD4	BV786	FITC/BB515/AF488/GFP	FITC/BB515/AF488/GFP
CD8	R718	PE/RYS86	PE/RYS86
CD45RA	PE-Cy7	APC/AF647	APC/AF647
CD197	BV711	BUV395	BUV395

Examples of fluorochromes that can be added as drop-ins using equivalent instrument and filter configurations


BUV: BD Horizon Brilliant™ Ultraviolet; BB: BD Horizon Brilliant™ Blue; R: BD Horizon™ Red; BV: BD Horizon™ Violet; AF: Alexa Fluor™

## Strategically designed to enable marker addition with minimal panel design and risk of resolution loss


### Optimal resolution of T cell subsets




### The drop-ins do not impact the backbone




### The backbone does not impact resolution of the drop-ins



### The recommended drop-ins do not impact each other



## Dissection of CD4+ Th cell subsets



## Methods

### Sample handling and staining

- Normal donor samples included lysed whole blood, freshly prepared PBMCs, and cultured and activated PBMCs.
- BD Horizon™ Brilliant Stain Buffer Plus is used at 10 µL/sample and included in all compensation controls, single-color controls and panels.
- Backbone panel reagents were used at 5 µL/test.
- Cells are pre-stained with CD197 (CCR7) BV711 reagent at 37 °C for 10 minutes. The same protocol was followed for the stain of any chemokine receptor in the CD4+ Th subset panel.
- Cells are stained with the full cocktail for 30 minutes.
- The BD Pharmingen™ Transcription Factor Buffer Set was used for the intracellular staining of FoxP3, according to manufacturer's instructions. Briefly, cells were stained with surface markers prior to fixation and permeabilization then intracellular markers are added for an additional 30-minute incubation at 4 °C.
- BD Pharmingen™ 7-AAD was added 10 minutes prior to acquisition, whereas BD Horizon™ Fixable Viability Stain (FVS) 620 was added during panel incubation in protein-free PBS.
- T cells were activated with Dynabeads™ Human T-Activator CD3/CD28 for 2 days.


### Instrument set up and acquisition

- BD FACSuite™ CS&T Research Beads were used to check QC of instrument performance
- Instruments were set up with optimal voltage application settings or lyse/wash settings for BD FACSLyric™ Cell Analyzer settings. Lyse/wash settings were adjusted to keep PE-Cy7 signal on scale
- Compensation was calculated with single-color control cells.


### Data Analysis

- FlowJo™ 10.8.1 Software was used for data analysis and Total Spread Matrix (TSM) table generation.
- Gates were drawn based on FMO controls.
- GraphPad Prism was used for graph generation.


## Identification of FoxP3+ Treg subsets



## Upregulation of inhibitory receptors



## Comprehensive T cell differentiation analysis



## Conclusions

- The human T cell backbone panel meets the criteria for a truly flexible and easy-to-expand flow cytometry panel:
  - ✓ Good resolution of major human T cell subsets
  - ✓ The backbone does not impact resolution of the recommended drop-in fluorochromes
  - ✓ The recommended drop-in fluorochromes do not impact the backbone resolution
  - ✓ The recommended drop-in fluorochromes do not impact each other's resolution
- The human T cell backbone panel is compatible with:
  - ✓ Several human sample types: lysed whole blood and fresh, cultured and/or activated PBMCs
  - ✓ Surface and intracellular stain protocols
  - ✓ DNA-binding dyes and Fixable Viability Stain for dead cell exclusion
  - ✓ Use with fluorescent protein GFP and streptavidin-PE, commonly used for detection of antigen-specific T cells or chimeric antigen receptors (CARs)
- The human T cell backbone panel was tested on different flow cytometers with different configurations (3, 4 and 5 lasers) and demonstrated consistent intra- and inter-instrument performance
- Up to five markers plus a viability dye can be added with minimal panel design effort without impact to population resolution and quantification, depending on instrument configuration
  - ✓ No redesign of the core backbone required in order to add new markers
  - ✓ No concerns about spillover, spread, compensation and co-expression when using the recommended drop-ins
  - ✓ Fluorochromes with appropriate brightness still need to be paired with markers based on antigen expression (bright-low, dim-high)
  - ✓ Fluorochrome with different brightness can be chosen for a given detector (e.g., dim FITC or bright BD Horizon Brilliant™ Blue 515 (BB515), dim V450 or bright BV421) thus providing further flexibility in panel design
- The human T cell backbone panel is strategically and prospectively designed to simplify the transition from five up to eleven color flow cytometry panels, leading to increased efficiency and biological insight
  - ✓ Double the number of markers analyzed in a single tube
  - ✓ Dive deeper into T cell biology through more comprehensive immunophenotypic and functional analyses