

Introduction

- Flow cytometry laboratories are increasingly integrating (semi) automated technologies into their routine workflows to meet the challenges of increasing workloads, assay complexity, limited human resources and regulatory requirements.
- The key areas in which automation can have the greatest impact are:
 - Reduction of Hands-On Time (HOT).
 - Reduction in Error Prone Tasks (EPT).
 - Process standardization and reproducibility.
- The BD FACSDuet™ Premium Sample Preparation System is a fully automated sample preparation device which integrates with the BD FACSLyric™ Flow Cytometer providing end to end sample processing including on-board sample washing and centrifugation.

Aims


- Compare fully automated sample processing with manually prewashed specimens completed on automation using:
 - Total sample processing time for both batched specimens and single specimens run consecutively.
 - Total hands-on time required during sample processing.
 - Number of error prone tasks.

Conclusions

- BD FACSDuet™ Premium Sample Preparation System provides reproducible and predictable sample processing times.**
- Hands-on time and error prone tasks are reduced when processes are fully (physically and digitally) automated.**
- Sample throughput can be increased by batched sample processing.**
- Maintaining some manual components, such as off-board pre-wash, may increase throughput but could increase variability, impact quality, and increase error prone tasks.**

Study Design and Method

Study Design



Specimen processing was completed at two sites using (a) full automation (BD FACSDuet™ Premium Sample Preparation System), inclusive of sample pre-washes, compared to (b) manual prewashing with process completed on automation (BD FACSDuet™ Premium Sample Preparation System). All manual prewash steps were matched on the automation. Total Processing Times (TPT), Hands-On Time (HOT), and number of Error Prone Tasks (EPT) were assessed on the following: (1) single specimen with two secondary tubes; (2) three consecutive worklists with a single specimen of two secondary tubes; and (3) batch of 8 specimens with 16 secondary tubes.

BD Biosciences provided materials and instruments for this study.

Disclaimers:

BD FACSDuet™ Premium Sample Preparation System and BD Flow Cytometers are Class I Laser Products

BD FACSDuet™ Premium Sample Preparation System is for In Vitro Diagnostic Use. Sample preparation for user-defined protocols and cocktail functions are for Research use Only, not for use in diagnostic or therapeutic procedures. BD FACSLyric™ Flow Cytometer is for Research Use Only with BD FACSuite™


BD, the BD Logo, BD FACSDuet, BD FACSLyric and BD FACSuite are trademarks of Becton, Dickinson and Company or its affiliates.

© 2023 BD. All rights reserved.

Analysis & Results

- Hands-On Time (HOT) and Error Prone Tasks are greatly reduced using the BD FACSDuet™ Premium Sample Preparation System for fully automated sample processing**
 - HOT was significantly lower (1.5% to 5.6% of TPT) compared to samples prewashed prior to completion on automation (4.6% to 12.7% of TPT).
 - Assuming the use of an integrated laboratory information system, fully automated sample processing reduces the number of tasks by 64.5% and reduces the Error Prone Tasks by 50%
 - Error Prone Tasks could be further reduced with the use of barcoded reagents on the BD FACSDuet™ Premium Sample Preparation System.


Condition tested	Tasks (Site1)	Error Prone (Site1)	%ERR
1 Sample OFF-Board prewash	31	4	12.9%
1 Sample ON-Board prewash	11	2	18.2%
Reduction with ON-Board prewash	64.5%	50.0%	



The total number of tasks and error prone tasks was determined for both sample handling procedures.


8 Specimen Worklists – 16 tubes; Samples were prepared in a single worklist; one worklist with pre-wash OFF-Board and one worklist with prewash ON-Board

Condition tested	TPT (Site1)	HOT (Site1)	%HOT (Site1)	TPT (Site2)	HOT (Site2)	TPT Difference - between sites	Site1 OFF-board: Increase throughput
8 Samples OFF-board prewash	2:17:09	0:17:26	12.7%	2:49:26	0:56:53	23.5%	
8 Samples ON-board prewash	3:15:03	0:10:58	5.8%	3:10:49	0:10:45	-2.2%	42.2%




Single Specimen Worklists – 2 tubes each; Samples were prepared in three consecutive worklists; one condition with pre-wash OFF-Board and one with prewash ON-Board

Condition tested	TPT (Site1)	HOT (Site1)	%HOT (Site1)	TPT (Site2)	HOT (Site2)	TPT Difference - between sites	Site1 OFF-board: Increase throughput
Sample 1 OFF-board prewash (Worklist 1)	1:59:47	0:07:56	8.2%	1:48:07	0:17:16	9.8%	
Sample 2 OFF-board prewash (Worklist 2)	2:36:17	0:07:11	4.6%	2:39:21	0:26:05	2.0%	
Sample 3 OFF-board prewash (Worklist 3)	2:31:23	0:09:21	6.2%	2:37:55	0:15:29	4.3%	
Sample 1 ON-Board prewash (Worklist 1)	1:55:34	0:04:58	4.3%	1:54:49	0:04:29	-0.6%	34.8%
Sample 2 ON-Board prewash (Worklist 2)	3:20:25	0:07:56	4.0%	3:20:01	0:04:19	-0.2%	
Sample 3 ON-Board prewash (Worklist 3)	3:23:45	0:02:59	1.5%	3:31:17	0:01:41	3.7%	



Legend



Data collection

The Lean component of this study used timers, paper logging, and video to capture the Total Process Time (TPT), Hands-On Time (HOT), and Error Prone Tasks (EPT) for time and motion. The time captured is from “Start of sample prep” to “ready for acquisition”.

- Using calibrated timers, video equipment was aligned with the instrumentation to ensure accuracy of record times (hh:mm:ss) for each step in the process to capture Total Process Time and Hands-On Time. Steps were also evaluated as to whether they were considered error prone.
- No patient identification was captured in documentation or by video equipment.
- Along with video equipment for tracking processes, paper records were made during the process in conjunction with the sites SOP.
- Lean specialists with a background in flow cytometry are crucial in identifying all steps and in the determination of error prone steps or deviations from SOP's that may lead to bias in the results.
- Laboratory staff performed tasks uninterrupted by the lean specialist to ensure there were no disruptions in the times observed or distractions from the SOPs.