

Antigen-specific T cell characterization by simultaneous detection of activation induced markers and intracellular cytokines: a multi-site testing experience

Gianluca Rotta¹, Valentina Achille², Melanie Duenas³, Marco Fernandez⁴, Federico Fondelli⁴, Arianna Gatti⁵, Sergio Haro Giron⁶, Martinez Rio Irene³, Diego J. Jiménez⁶, Jorge Monserrat⁶, Verena Pfeifer⁷, Josefa Pichler⁸, Barbara Prietl⁷, Harald Sourij⁹, Chiara Uras¹⁰, Sara Vlah¹¹, Daniela Fenoglio¹²

¹ BD Biosciences Europe, ² Cytogenetics and Molecular Genetics Laboratory ASST Ovest Milanese - Legnano Hospital - Italy, ³ BD Biosciences - Spain, ⁴ Flow Cytometry Core Facility Germans Trias i Pujol Research Institute - Spain, ⁵ Transfusion Center ASST Ovest Milanese - Legnano Hospital - Italy, ⁶ Department of Medicine and Medical Specialties, University of Alcalá, Alcalá de Henares - Spain, ⁷ Division of Endocrinology and Diabetology, Medical University of Graz, Center for Biomarker Research in Medicine GmbH, Graz - Austria, ⁸ BD Biosciences Austria, ⁹ Division of Endocrinology and Diabetology, Medical University of Graz, ¹⁰ Department of Experimental Medicine, Centre of Excellence for Biomedical Research, University of Genoa - Italy, ¹¹ BD Biosciences Italy, ¹² Department of Internal Medicine, University of Genoa, Biotherapy Unit, IRCCS Ospedale Policlinico San Martino - Genoa, Italy

Abstract

Introduction. During the COVID-19 pandemic, the capacity of several SARS-CoV2 variants to escape the vaccine-induced neutralizing antibody response, together with the observation that T cell mediated immunity was essentially conserved, led to an increased interest in the characterization of antigen-specific T cells by flow cytometry. In this context, we recently published an optimized protocol for simultaneous detection of T cell activation-induced markers and intracellular cytokines, which is based on overnight PBMC stimulation with Spike peptide pool and subsequent flow cytometry characterization using a 12-marker panel (Live-Dead, CD3, CD4, CD8, CD45RA, CCR7, CD137, CD69, CXCR5, IFN γ , TNF α , IL-2)(see references). With the aim of understanding if this protocol could be successfully reproduced in different laboratories, we established a network of five testing sites who characterized a total of 40 SARS-CoV2 reactive subjects accounting for 200 acquired samples.

Methods. Each site worked with a different flow cytometry platform: four analyzers and one cell sorter, three out of five instruments being equipped with a yellow-green laser. Operators were trained and supported during the first experiment execution, a common acquisition template was distributed, instrument settings were harmonized using hard-dyed beads, and data analysis was centralized (Figure 1). Five specific aspects were investigated: 1) the capacity to successfully report data for each analyzed sample; 2) the impact on Stimulation Index (S.I. = ratio between activated T cells in stimulated and unstimulated samples) of intracellular (ic) vs surface (s) CD4+ and CD8+ staining; 3) the impact on S.I. of aCD28 co-stimulation during PBMC activation; 4) the consistency of analysis when different operators evaluated the same raw data files; 5) the estimated limit of detection of the assay.


Results. Two donors out of 40 were discarded for poor sample quality and three were not reported after analysis due to low cell counts; in addition, three samples required a modification in the gating strategy to remove unwanted background. CD4+ and CD8+ intracellular vs surface staining were assessed on nine donors, revealing a lower ic-CD4+ vs s-CD4+ S.I. (ic-CD4+ S.I. = 5.3; s-CD4+ S.I. = 8.7; p-Value = 0.036) (Figure 2). No significant differences were detected for the CD8+ subset. Interestingly, samples activated without aCD28 (n=26) revealed a significantly higher CD4+ S.I. (9.4 vs 4.3; p-Value = 0.001) and cytokine-producing CD4+ S.I. (13.6 vs 6.4; p-Value = 0.002) (Figure 3), which was driven by a 54% reduction of background activation in unstimulated samples for the total activated CD4+ subset and by a 43% reduction of background activation for the cytokine-producing CD4+ subset (Figure 4). Corresponding stimulated samples were also characterized by lower activation yet to a milder extent (18% and 16%, respectively) (Not shown). No relevant variations were reported for the CD8+ compartment.

When different operators analyzed the same raw data from three stimulated samples, coefficients of variation for reported subsets were inversely related to the abundance of cells. CD4+ and CD8+ being respectively 65.8% and 22.6% of CD3+ were reported with %CV of 2.6 and 9.4, whereas activated CD4+ and CD8+ being respectively 0.11% and 0.30% of parental T helper and cytotoxic cells were reported with %CV of 18.6 and 30.5.

10⁶ PBMCs were stimulated per sample, leading to an average acquisition of 21,611 and 11,720 live CD4+ and CD8+ cells/sample. Considering the requirement of 20 recorded events to detect a given subset, the average lower limit of measurable activation was 0.09% for CD4+ and 0.17% for CD8+. One stimulated sample out of 26 displayed CD4+ activation below the limit and two stimulated samples displayed CD8+ activation below the limit. S.I. \geq 3 was required to define responder samples: CD4+ response was detected in 20 out of 25 samples and CD8+ response was detected in 8 out of 24 samples.

Conclusions: The method was successfully implemented in all the testing sites, use of s-CD4 and s-CD8 staining as well as avoidance of co-stimulation with aCD28 are recommended. Limit of detection is estimated at 0.09% for CD4+ and 0.17% for CD8+ events.

Figure 1 - Characterization of different donors performed by diverse operators on diverse instruments provides similar staining patterns



A Flow cytometry panel overview


B Instrument setting harmonization

Instrument settings were harmonized by means of BD FACSDiva™ CS&T Research Beads. The same fluorescent target values were set in all instruments for every channel and beads were acquired before every working session. The panel shows channel specific MFIs and %CVs across instruments and acquisitions.

C Exemplary CD4+ T cell staining pattern

The panel shows exemplary CD4+ T cell staining patterns from independent SARS-CoV-2 reactive samples (frozen PBMCs) activated with Spike peptide pools. Cells were processed by different operators working in different sites and acquired with diverse platforms. Each operator was trained on the protocol and the analysis was centralized. Staining patterns always looked similar during the testing phase and could be analyzed using the same gating strategy pending small sample specific gate adjustment.


Figure 2 - CD4 and CD8 markers do not require intracellular (ic) staining



CD4+ and CD8+ T cell Stimulation Index with surface and intracellular staining

9 samples were stimulated across the different sites with Spike peptide pools and subsequently stained with the 12-color panel using surface or intracellular CD4 and CD8 antibodies. Stimulation Index (S.I.) was calculated as the ratio between Activation Induced Marker positive (AIM+) cells in stimulated and unstimulated samples. The use of intracellular CD4 and CD8 markers did not result in higher Stimulation Index, rather ic-CD4 staining resulted in slightly higher background activation thus lowering the related Stimulation Index.


Figure 3 - CD4+ Stimulation Index is higher without using aCD28 co-stimulation



CD4+ Stimulation index with or without aCD28 co-stimulation

26 samples were stimulated across the different sites with Spike peptide pools in the presence or absence of aCD28 used as co-stimulus. Unstimulated samples were treated accordingly with or without aCD28. Multicolor flow cytometry revealed a significantly higher Stimulation Index when samples were not co-stimulated with aCD28. The pattern was conserved for total activated and for cytokine producing CD4+ T cells.

Figure 4 - Absence of aCD28 results in lower background activation of unstimulated samples



CD4+ background activation

CD4+ T cell background activation was lower in the absence of aCD28 co-stimulation, either for total activated CD4+ T cells (AIM+) and for the cytokine-producing subset (Cyto+ AIM+).

Methods

10⁶ cells/sample were seeded per well in a 96-u bottom plate and stimulated with 1 μ g/ml of PepTivator® SARS-CoV-2 Prot S, Prot S1, Prot S+ or with 1 μ g/ml of PepTivator® of CMV pp65 and CMV IE-1 (Miltentyi). Monoclonal antibodies anti-CD28 at 1 μ g/ml were added as co-stimulus. After 4h of incubation at 37°C, 5% CO₂, the protein transport inhibitor BD Golgi-Plug™ (Brefeldin A, 1 μ l/ml) was added for a further incubation of 20 hours. Detailed methods are available at

<https://www.sciencedirect.com/science/article/pii/S002217592300025X?via%3Dihub#s0055>

Figure 5 - Lower limit of detection and reproducibility of analysis

Lower limits of detection are
0.09% for AIM+ CD4+
and 0.17% for AIM+ CD8+

10⁶ PBMCs were stimulated per sample, leading to an average acquisition of

- CD4+ average acquired events = 21,611
- CD8+ average acquired events = 11,720.

Requirement of 20 recorded events to detect activated cluster was considered according to F. Buccisano et al. Haematologica. 2022 Dec 1; 107(12): 2823-2833. doi: 10.3324/haematol.2021.279777

Different operators provided consistent results when analyzing the same raw FCS data (N=3)

%CV for reported CD4+ = 2.6
%CV for reported CD8+ = 9.4

%CV for reported AIM+ CD4+ = 18.6
being AIM+ CD4+ = 0.11% of total CD4+

%CV for reported AIM+ CD8+ = 30.5
being AIM+ CD8+ = 0.3% of total CD8+

Conclusions

This method was adopted in five different laboratories and the 12-color panel provided meaningful results on 38 samples out of 40 when acquired on different flow cytometry platforms. The use of intracellular staining for CD4 and CD8 markers is not required and the co-stimulation with aCD28 increased T cell background activation thus reducing the final CD4+ Stimulation Index.

Class 1 Laser Product. For Research Use Only. Not for use in diagnostic or therapeutic procedures. BD, the BD logo and FACS are the trademarks of Becton, Dickinson or its affiliates. © 2024 BD. All rights reserved. BD-120757(v1.0) 0524