

# BD FACSPresto™ Near-Patient CD4 Counter

Complete CD4 testing solution for advanced HIV disease

## Simple workflow with high precision and accuracy



Choice of finger stick or venipuncture draw



18 minute incubation occurs outside device at room temperature



4 minute read

Quick results in 22 minutes total time

- Absolute CD4, %CD4, and Hb in a single test
- Detects wide range of CD4 counts (50–4,000 cells/ $\mu$ L blood)
- Results on screen, in print or exported using USB

## Easy to use for resource limited settings



Easy to use, small and portable



Onboard videos and timers, language independent



No cold chain required with long shelf life



Built-in battery to ensure operation during power interruptions



Built-in quality controls, no need for additional costly external controls

## Affordable



Accessible pricing  
Compared to systems currently on the market



No added cost  
Instrument placement with reagent rental agreement (minimum 3 years, 1,500 tests per year)



Included benefits

- In country training and on-site support
- Swap and repair offerings
- <72 hours guaranteed response

For In Vitro Diagnostic Use.

For E.U. under IVDD and IVDR, the BD FACSPresto™ Near-Patient CD4 Counter and BD FACSPresto™ Cartridges are for In Vitro Diagnostic Use.

The BD FACSPresto™ Cartridge is an in vitro diagnostic device bearing a CE marking and is CE certified by BSI Group the Netherlands B.V. (Notified Body Number : 2797).

BD Life Sciences, San Jose, CA 95131, USA

[bdbiosciences.com](http://bdbiosciences.com)

BD, the BD Logo and BD FACSPresto are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.  
BD-65713 (v1.0) 0922 / Formerly 23-22135-00

