



# BD FACSVia™ System

Product list



# BD FACSVia™ System

## Instrument

Description	Cat. No.
BD FACSVia™ Flow Cytometer (includes the accessory box)	656874

## Hardware

Description	Cat. No.
BD FACSVia™ PC Bundle (includes BD FACSVia™ clinical software)	660508 US/660509 EU
23-inch NEC Multisync Monitor	652792
BD FACSVia™ Clinical Software	660562

## Options

Description	Cat. No.
BD FACSVia™ Loader (for sample automated loading) (includes two 24-tube racks with key)	656875
BD Hand-Held Barcode Reader (ISBT 128 supported)	657929

## Options for LIS connectivity

Description	Cat. No.
BD FACSLink™ Next Generation Base Software for BD FACSVia™ System	660733
BD FACSLink™ Computer Assembly	660770
BD FACSLink™ Next Generation Connection License for BD FACSVia™ System (for optional multiple instruments)	660735

## Option BD FACSVia research software

Description	Cat. No.
BD FACSVia™ Research Software (for user defined protocols)*	660563

## Fluids and beads

Description	Cat. No.
BD™ Sheath Additive (10 bottles of 5 mL per box) Add one bottle (5 mL) of BD Sheath Additive to 1 liter of 0.2-µm filtered DI water.	660584
BD™ Detergent Solution Concentrate (bottle of 15 mL) Add 3 mL of BD Detergent Solution Concentrate to 197 mL of filtered DI water. Use the working solution within 2 weeks.	660585
BD™ FACSClean Solution (cubitaiver of 5 L)	340345
BD™ Extended Flow Cell Clean Solution	660586
BD™ CS&T Beads 50 or 150 test (IVD, specifically for BD FACSVia).	656504 (50 test) CEIVD 656505 (150 test) CEIVD

## Customer maintenance

Description	Cat. No.
2-Month Preventative Maintenance Kit	660634
12-Month Preventative Maintenance Kit	661331

## Consumables and spare parts

Description	Cat. No.
Filters for BD™ Sheath Additive, BD™ Detergent Solution Concentrate and BD™ FACSClean Bottles (3 pieces)	660322
Peristaltic Pump Tubing (2 pieces)	659522
Waste Bottle - 2 L (includes cap and sensor)	659605
BD™ Sheath Additive Bottle - 2 L (includes cap and sensor)	661277
BD™ FACSClean Bottle 250 mL (includes cap and sensor)	661275
BD™ Detergent Solution Concentrate Bottle 250 mL (includes cap and sensor)	661276

\*For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Class 1 Laser Product.

For In Vitro Diagnostic Use. CE-IVD. Not available in USA.

23-17579-00

BD Life Sciences, San Jose, CA, 95131, USA

[bdbiosciences.com](http://bdbiosciences.com)

