

---

BD Life Sciences

# BD Accuri™ Product List

---



## BD Accuri™ Product List

### Systems and Software

DESCRIPTION	NOTES	CAT. NO.
BD Accuri™ C6 Flow Cytometer System, Software, and Workstation Computer	Includes standard optical filters, reagent containers, complimentary BD Accuri™ C6 flow cytometer starter kit, and one-year warranty.	653118
BD Accuri C6 Flow Cytometer System, BD CSampler™, Software, and Workstation Computer	Includes standard optical filters, reagent containers, BD Accuri C6 flow cytometer starter kit, and one-year warranty for BD Accuri C6 and BD CSampler.	653119
BD Accuri™ C5 Flow Cytometer System, Software, and Workstation Computer	Includes 488-nm laser, 2 light scatter and 3 fluorescence detectors, complimentary BD Accuri flow cytometry starter kit, and one-year warranty.	657214
BD CSampler: 96-well plate and 24-tube rack Automated Sampler for BD Accuri C6 Flow Cytometer	Includes 24-tube rack, Enhanced Analysis Module, and one-year warranty.	653124

### Optional System Modules

DESCRIPTION	NOTES	CAT. NO.
Selectable Lasers Module for BD Accuri C6 Flow Cytometer	Facilitates standard three detectors off blue laser, and one off red laser. Allows for four blue detectors, or two blue and two red detectors. Includes Activation Key and Optical filters: 780/60, 610/20, 630/30.	653126
User Tracking Module for BD Accuri C5 and C6 Flow Cytometers	For the creation of unique user names and passwords to record the duration each operator uses the cytometer. A log file is created, allowing the preparation of billing reports for multi-user or core laboratories.	653127

### Flow Cytometry Analysis Software

DESCRIPTION	NOTES	CAT. NO.
BD Accuri C6 Analysis Software for PC or Mac®	For use with data exported from the BD Accuri C6	653122
BD CSampler Analysis Software for PC or Mac	For use with data exported from the BD Accuri C6 with BD CSampler	653123
FCS Express™ 4 Flow Cytometry RUO Professional, Single License	With C6 file importer	653129
MultiCycle AV DNA Cell Cycle Analysis Software, Single License	Requires a valid FCS Express v4 Professional license	653130

### Extended Warranties

DESCRIPTION	NOTES	CAT. NO.
Service Contract for BD Accuri C6 Flow Cytometer (per year)	Includes complimentary 1-year Maintenance Kit (653149)	660761
Service Contract for BD CSampler (per year)		660763
Service Contract for BD Accuri C5 Flow Cytometer (per year)	Includes complimentary 1-year Maintenance Kit (653149)	660762
Service Contract for BD Accuri C6 SORP Flow Cytometer (per year)	Includes complimentary 1-year Maintenance Kit (653149)	660768

## BD Accuri Product List

### Maintenance, Quality Assurance, and Service

DESCRIPTION	NOTES	CAT. NO.
Spherotech 8-Peak Validation Beads (FL1–FL3)	4.0 mL	653144
Spherotech 6-Peak Validation Beads (FL4)	2 x 4.0 mL	653145
Peristaltic Pump Tubing (Set of 2)		653146
Sheath, Decon and Cleaner Bottle Filters (1 each)		653147
In-line Sheath Filter (1)		653148
BD Accuri Flow Cytometer 1-year Maintenance Kit	Contains (12) Peristaltic Pump Tubing, (6) Sheath Bottle Filters, (12) Decon/Cleaner Bottle Filters, and (6) In-Line Sheath Filters	653149
BD Accuri Flow Cytometer 2-month Maintenance Kit	Contains (2) Peristaltic Pump Tubing, (1) Sheath Bottle Filter, (2) Decon/Cleaner Bottle Filters, and (1) In-Line Sheath Filter	653150
Installation and New User Training Service	The majority of BD Accuri C6 users have found it easy to install the BD Accuri C6 flow cytometer using the provided guidance documents and are typically up and running in 2 hours or less. If desired, we can provide onsite, professional instrument set up and installation with 4 hours of basic instrument training. Price includes BD Accuri employee travel.	653151
New User Flow Cytometry Training Service	Service to include 4–6 hours of training provided onsite (includes travel)	653283

### Reagents and Kits

DESCRIPTION	NOTES	CAT. NO.
Decontamination Concentrate Solution	(Makes 200 mL)	653154
Decontamination Concentrate Solution, Package of 5	(Makes 1.0 L)	653155
Bacteriostatic Concentrate Solution (for sheath fluid)	10 vials/package (treats a total of 10 L of sheath)	653156
Cleaning Concentrate Solution	(Makes 1.0 L)	653157
BD Accuri Flow Cytometer Fluid Kit	BD Accuri flow cytometer fluid kit contains concentrate for making 1 liter each of Cleaning and Decontamination Solutions and Bacteriostatic Concentrate for treating 10 L of sheath, as well as 8 mL of our Extended Flow Cell Clean Solution.	653158
Extended Flow Cell Clean Solution (8 mL)		653159
BD Accuri Flow Cytometer Starter Kit	BD Accuri flow cytometer starter kit contains concentrate for making 1 liter each of Cleaning and Decontamination Solutions and Bacteriostatic Concentrate for treating 10 L of sheath. It also contains 4 mL of each of the 6- and 8-peak Spherotech Validation Beads, as well as 8 mL of our Extended Flow Cell Clean Solution.	653160

### Replacement Parts and Accessories

DESCRIPTION	NOTES	CAT. NO.
24-Tube Rack for BD CSampler		653161
BD CSampler Protector Pelican Case	Includes packaging foam	653162
BD Accuri Flow Cytometer Protector Pelican Case	Includes packaging foam	653163
Waste Bottle (Serial Numbers <2400)*		653164
Sheath Bottle (Serial Numbers <2400)*		653165
Cleaning Solution Bottle (Serial Numbers <2400)*		653166
Decontamination Bottle (Serial Numbers <2400)*		653167
2-L Waste Bottle (Serial Numbers ≥2400)*		653168
2-L Sheath Bottle (Serial Numbers ≥2400)*		653169
Cleaning Solution Bottle (Serial Numbers ≥2400)*		653170
Decontamination Bottle (Serial Numbers ≥2400)*		653171

\*Products include complete level sensor assembly (bottle, cap, and level sensor).

## BD Accuri Product List

### Optical Filters

DESCRIPTION	CAT. NO.
FL1 99% Attenuation Filter	653172
FL1 90% Attenuation Filter	653173
FL2 99% Attenuation Filter	653174
FL2 90% Attenuation Filter	653175
FL3 99% Attenuation Filter	653176
FL3 90% Attenuation Filter	653177
FL4 99% Attenuation Filter	653178
FL4 90% Attenuation Filter	653179
Standard Optical Filter 533/30 nm (FL1)	653180
Standard Optical Filter 585/40 nm (FL2)	653181
Standard Optical Filter >670 nm LP (FL3)	653182
Standard Optical Filter 675/25 nm (FL4)	653183
Optical Filter 510/15 nm	653184
Optical Filter 565/20 nm	653185
Optical Filter 610/20 nm	653186
Optical Filter 780/60 nm	653187
Optical Filter 540/20 nm	653528

Please note flow cytometer serial number prior to contacting Customer Service to ensure purchase of correct optical filters.

Applicable Terms and Conditions may be found on the BD Accuri Cytometers website at [bdbiosciences.com](http://bdbiosciences.com).

Class 1 Laser Product.

Trademarks are the property of their respective owners.

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



BD Life Sciences  
Biosciences  
[bdbiosciences.com](http://bdbiosciences.com)